



Glean Urodynamics System

Owner's Manual

Operation, Care, and Maintenance



Bright Uro, Inc.
3 Goddard
Irvine, CA 92618 USA
Tel: +1 (949) 216-0873
Fax: +1 (949) 998-4731
Email: info@brighturo.com



Bright Uro, Inc.
3 Goddard
Irvine, CA 92618 USA
Tel: +1.949.216.0873
Fax: +1 (949) 998-4731
Email: info@brighturo.com

www.brighturo.com

Glean Urodynamics System and GUS are trademarks of Bright Uro, Inc. Bluetooth is a registered trademark of Bluetooth SIG. Windows is a registered trademark of Microsoft Corp. Other names may be trademarks of their respective owners.

ALL RIGHTS RESERVED. IMAGES/RENDERINGS OF BRIGHT URO, INC. DEVICES IN THIS PUBLICATION ARE PROPERTY OF BRIGHT URO AND NO PART OF THIS PUBLICATION MAY NOT BE REPRODUCED WITHOUT THE PRIOR WRITTEN PERMISSION OF BRIGHT URO.

© Copyright 2023–20XX Bright Uro, Inc.

Document Number: LBL-XXXX Revision:

0

Release Date: 202X-XX



Bright Uro, Inc.
3 Goddard
Irvine, CA 92618 USA

R_x ONLY

GLEAN URODYNAMICS SYSTEM	I
ABOUT THIS MANUAL	7
SYMBOLS	7
1 INTRODUCTION.....	8
1.1 DEVICE DESCRIPTION	8
1.1.1 Sensor.....	8
1.1.2 Insertion Tool.....	9
1.1.3 Uroflowmeter.....	9
1.1.4 Software – Glean Mobile Apps (Clinician and Patient) and Glean Web App.....	9
1.2 GETTING STARTED	10
1.2.1 How Supplied.....	10
1.2.2 Device Inspection.....	11
1.2.3 Device Storage	11
1.3 LEARNING ABOUT THE GLEAN URODYNAMICS SYSTEM (GUS)	11
1.4 INTENDED USE / INDICATIONS FOR USE.....	11
1.5 CONTRAINDICATIONS.....	11
1.6 TARGET USERS	11
1.7 WARNINGS AND PRECAUTIONS.....	12
2 ADVERSE EVENTS/RESIDUAL RISKS	13
2.1 POTENTIAL ADVERSE EVENTS/RESIDUAL RISKS.....	13
2.2 DEVICE-RELATED ADVERSE EVENTS REPORTING	14
3 PATIENT COUNSELING INFORMATION	14
4 GLEAN URODYNAMICS SYSTEM—CARE AND MAINTENANCE	14
4.1 GENERAL CARE, CLEANING, AND PREVENTIVE MAINTENANCE	14
4.1.1 Caring for the GUS Uroflowmeter.....	14
4.2 BATTERY—CHARGING AND PREVENTATIVE MAINTENANCE	14
4.2.1 Charging the Battery	14
4.4 TREATING AND DISPOSING OF PRODUCT AFTER USE	15
4.5 ENVIRONMENTAL CONSIDERATION OF WASTE DISPOSAL	15
4.6 PREVENTIVE MAINTENANCE—CHECKING CALIBRATION	15
5 SET UP THE GLEAN URODYNAMICS SYSTEM (GUS)	15
6 EQUIPMENT STATUS INDICATORS AND BUTTONS	16
6.1 LED LIGHTS	16
6.1.1 Sensor LED Patterns	16
6.1.2 Uroflowmeter LED Patterns	16
7 SOFTWARE FEATURES AND FUNCTIONS.....	16
7.1 ACCESSING THE GLEAN APPS	16
7.1.1 Access the Glean Web App	17
7.1.2 Install the Glean Mobile Apps	17
7.2 LOGIN TO THE GLEAN MOBILE OR WEB APPS	17
7.3 CREATE USER ACCOUNTS	18
7.4 DELETE USER ACCOUNTS.....	18
7.5 RESET PASSWORD	20
7.5.1 Admin	20
7.5.2 User.....	21
7.6 ADD A PATIENT PROFILE IN THE GLEAN ADMIN PORTAL	21

8 HOW TO RUN TESTS – CMG/PF TEST, UROFLOW TEST, AND DATA ANALYSIS	23
8.1 CMG/PF TEST.....	23
8.1.1 Prepare the patient for aseptic insertion.....	23
8.1.2 Instill lubricant in the urethra.....	23
8.1.3 Prepare the Sensor for data collection.....	23
8.1.4 Load the Sensor in the Sheath.....	28
8.1.5 Deploy the Sensor in the bladder.....	28
8.1.6 Log events using the Glean Mobile App (Clinician).....	29
8.1.7 Instruct the patient to log symptoms during ambulatory monitoring.....	33
8.1.8 Train the Patient to use the Glean Mobile App (Patient).....	34
8.1.9 Run a Uroflow Test.....	36
8.1.10 Measure PVR using preferred method.....	43
8.1.11 Prepare the patient for Sensor removal.....	43
8.1.12 Remove the Sensor from the bladder.....	43
8.1.13 Urethral Pressure Profile (UPP).....	43
8.1.14 Download data from the Glean sensor.....	43
8.2 ANALYZE THE DATA USING THE GLEAN WEB APP	47
8.2.1 Login to the Glean Web App (Clinician).....	47
8.2.2 Select the desired patient and study.....	47
8.2.3 Analyze Urodynamics data.....	48
8.2.4 Adjust the x-axis for Urodynamics data.....	48
8.2.5 Adjust the y-axis for Urodynamics data.....	50
8.2.6 Select multiple events of similar type.....	50
8.2.7 Draft notes for interpretation, assessment, and plan.....	50
8.2.8 Export the Urodynamics report.....	51
8.2.9 Utilize filtering to adjust view of Urodynamics data.....	53
8.2.10 View the Glean patient history.....	53
8.2.11 View detailed Urodynamic study data.....	54
8.2.12 Adjust Glean Web App settings.....	55
9 CALIBRATION	56
10 FAQ	57
11 TROUBLESHOOT	57
12 APPENDICES.....	60
APPENDIX A: TECHNICAL DATA.....	60
APPENDIX B: CLASSIFICATIONS	61
APPENDIX C: SYMBOLS AND LABELING.....	62
APPENDIX D: ELECTROMAGNETIC COMPATIBILITY (EMC)	63
APPENDIX E: END-USER SOFTWARE LICENSE AGREEMENT	69
APPENDIX F: GLOSSARY	69
APPENDIX G: ACRONYMS	72
13 REFERENCES.....	74

LIST OF FIGURES

Figure 1. Sensor	8
Figure 2. 18 Fr. Insertion Tool - Advancer and Sheath (male-left; female-right)	9
Figure 3. Uroflowmeter	9
Figure 4: Software (Clinician App, Patient App, Web App)	10
Figure 5. Urine Collection Cup Placement	15
Figure 6. Glean Mobile App Bluetooth Connection	17
Figure 7. Glean Mobile App and Web App Login Page	17
Figure 8. Navigate to Account Page	18
Figure 9. Create an Account	18
Figure 10. Enter Account User Information	18
Figure 11. Navigate to Account Page	19
Figure 12. View User Account	19
Figure 13. Delete User Account	19
Figure 14. Confirm Deletion of User Account	19
Figure 15. Navigate to Account Page	20
Figure 16. View User Account	20
Figure 17. Reset Password	20
Figure 18. Reset Password Email Successfully Sent	21
Figure 19. Forgot Password	21
Figure 20. Send Email to Reset Password	21
Figure 21. Patient Account Page	22
Figure 22. Create Patient Account	22
Figure 23. Enter Patient Information	22
Figure 24. Start a New Study	23
Figure 25. Select a Patient	24
Figure 26. Confirm Patient Information	24
Figure 27. Start Study	24
Figure 28. Add Events or Devices	25
Figure 29. Add Sensor	25
Figure 30. Power On the Sensor	26
Figure 31. Scan the QR Code on the Pouch	26
Figure 32. Enter PIN and Setup Sensor	27
Figure 33. Calibrate the Sensor	28
Figure 34. Sensor Deployment Complete or Failure	29
Figure 35. Add Events or Devices	30
Figure 36. Add Events (Clinician)	30
Figure 37. Select Appropriate Event (Clinician)	31
Figure 38. Confirm Event Details (Clinician)	32
Figure 39. Event Timer	33
Figure 41. Edit Event (Clinician)	33
Figure 41. Add Event (Patient)	34
Figure 42. Select Event Type (Patient)	34
Figure 44. Enter Event Details (Patient)	35
Figure 45. Edit Event (Patient)	35
Figure 46. Delete Event (Patient)	36
Figure 46. Start a New Study	37
Figure 47. Select a Patient, Confirm Patient Information, and Start Study	37
Figure 48. Add Events or Devices	38
Figure 49. Add Uroflowmeter	38
Figure 50. Power on the Uroflowmeter	39
Figure 51. Scan QR Code on Uroflowmeter	39
Figure 53. Enter Uroflowmeter PIN	40
Figure 54. Associate Uroflowmeter with Patient	40
Figure 55. Start Uroflowmetry	41
Figure 55. Stop Uroflowmetry	41
Figure 57. Download Uroflowmetry Data	42
Figure 58. Confirm Uroflowmetry Data	43
Figure 59. Add Events or Devices	44
Figure 60. Select Download Sensor Data	44
Figure 61. Connect Sensor to Download Data	45

Figure 61. Download Sensor Data	46
Figure 62. Confirm Sensor Data.....	46
Figure 64. Studies Overview	47
Figure 64. Select Patient Studies	47
Figure 65. View Patient Studies	47
Figure 66. View Study Data.....	48
Figure 67. Drag Zoom Box.....	48
Figure 68. Adjust Sliding View Window.....	49
Figure 69. Select Event of Interest.....	49
Figure 70. View Event of Interest	49
Figure 71. Uroflowmetry and Voiding Diary Tabs	50
Figure 72. Select Events of Similar Type	50
Figure 73. Unselect Events of Similar Type	50
Figure 74. Draft Primary Interpretation.....	50
Figure 75. Draft Assessment and Plan	51
Figure 76. Copy Notes	51
Figure 77. Save Notes.....	51
Figure 78. Complete Urodynamics Report.....	52
Figure 79. View Urodynamics Report	52
Figure 80. Print Urodynamics Report.....	52
Figure 81. View Filters.....	53
Figure 83. Select Filters	53
Figure 83. Select Patient Record and Study	54
Figure 84. View Glean Patient History	54
Figure 85. Select Patient Record and Study	54
Figure 86. View Study Data.....	55
Figure 88. Glean Web App Settings.....	55
Figure 89. Adjust Glean Web App Settings.....	55
Figure 90. Confirm Settings Changes	56

LIST OF TABLES

Table 1. Required Equipment	10
Table 2. Device LED Light Locations	16
Table 3. Uroflowmeter Device State LED Patterns	16
Table 4. Uroflowmeter Battery State LED Patterns.....	16
Table 5. Troubleshooting.....	59
Table 6. Error Messages	60
Table 7. GUS Specifications	61
Table 8. GUS Classifications.....	61
Table 9. Symbols Glossary	63
Table 10. Label Location	63
Table 11. Electromagnetic Compatibility	64
Table 12. Table 1 Requirements—Electromagnetic Environment	65
Table 13. Table 2 Requirements—Electromagnetic Environment—Guidance.....	66
Table 14. Table 4 Requirements—Electronic Environment—Guidance	67
Table 15. Transmitter and Receiver Product Specifications	68
Table 16. Minimum Separation Distance	68

ABOUT THIS MANUAL

SYMBOLS

This manual provides important information to help in understanding the features and safe use of the device. The symbols outlined here highlight helpful tips and important cautions that will aid in guiding the reader through the manual.



Caution/warning symbol describes information that the user needs to know to prevent minor injury or product damage.



Important symbol describes important information about using the device.



Note symbol describes additional information about the device.

1 INTRODUCTION

1.1 DEVICE DESCRIPTION

The Glean Urodynamics System (GUS) is indicated for standard Urodynamic tests such as Uroflow (UF), Cystometrogram (CMG), Urethral Pressure Profile (UPP), and Micturition Studies (MS).

GUS consists of the following three physical component elements: Sensor, Insertion Tool, and Uroflowmeter, as well as the following three software applications: Glean Mobile App (Clinician), Glean Mobile App (Patient), and Glean Web App. The patient may use the Glean Mobile App as a digital voiding diary, logging fluid input, leakage, urgency, and other urological symptoms. The clinician may use the Glean Mobile App to prepare the Sensor for insertion, log symptoms, and download data. The Glean Web App may be used by clinicians to view and analyze data.

The Sensor can be inserted through the urethra into the bladder using the Insertion Tool. Once inserted, the Sensor has a Removal String that hangs out of the urethra to enable removal of the Sensor. The Sensor may stay in the bladder for the entire duration of monitoring while collecting data. The Sensor stores data that may be wirelessly transmitted to the Glean Mobile App (Clinician) once it is removed from the body.

The Uroflowmeter is used to measure voided volume and flow. The Glean Mobile App (Clinician) wirelessly receives data from the Uroflowmeter after the patient has completed a voiding cycle.

Data from the Glean Mobile App (Clinician and Patient) is synchronized wirelessly in the cloud and made available for review, analysis, and interpretation by trained clinicians using the Glean Web App. After the clinician has completed analysis and interpretation, they may use the Glean Web App to generate the report documenting the Urodynamic findings.

1.1.1 Sensor

- The Sensor (Figure 1) is a long, flexible tube with rounded ends and a removal string.
- The Sensor is designed to be inserted into the bladder through the urethra using the Insertion Tool.
- The Sensor contains a battery and flexible electronic circuit board with a microprocessor, software, pressure sensor and memory to store data.
- The Sensor is designed to curl into a circular shape once inserted in the bladder to ensure the Sensor stays in the bladder until removal is desired.
- When desired, the clinician may remove the Sensor by gently pulling on the Removal String. This will pull the Sensor out of the body through the urethra.
- The Sensor can collect data for the duration of Urodynamic monitoring.
- The Sensor is a one-time use disposable device and designed for use under the supervision of a trained clinician.



Figure 1. Sensor

1.1.2 Insertion Tool

- The Insertion Tool (Figure 2) is used to insert the Sensor in the patient's bladder.
- The Insertion Tool is comprised of two components: the Sheath and the Advancer.
- Once the Sensor is placed in the bladder, the Sheath and Advancer may be removed leaving the Removal String hanging out of the urethra.




Figure 2. 18 Fr. Insertion Tool - Advancer and Sheath (male-left; female-right)

1.1.3 Uroflowmeter

- The Uroflowmeter (Figure 3) measures the voided volume and flow of urine when a patient voids into the urine collection cup.
- The Uroflowmeter is designed to work with commonly available commodes and funnels that may assist users in properly collecting the urine.



Figure 3. Uroflowmeter

1.1.4 Software – Glean Mobile Apps (Clinician and Patient) and Glean Web App

- The Glean Mobile App (Clinician) is used by the clinician to prepare the Sensor for insertion.
- During the period of monitoring, the patient may use the Glean Mobile App (Patient) to complete a digital voiding diary. The patient may log fluid input, leakage, urgency, and other urological symptoms.
- After monitoring is complete, the clinician may use the Glean Mobile App (Clinician) to download data from the Sensor.
- The Glean Mobile App (Clinician) is used to download data from the Uroflowmeter.
- Data from the Glean Mobile App (Clinician and Patient) is synchronized in the cloud and available for review, analysis, and interpretation by a trained clinician using the Glean Web App.
- Once complete, the clinician may use the Glean Web App to generate a Urodynamics report.




Figure 4: Software (Clinician App, Patient App, Web App)

1.2 GETTING STARTED

To prepare the GUS before performing any Urodynamics procedures:

1. Remove the Uroflowmeter from the packaging and ensure it powers on by pressing and holding the Button LED on the front of the Uroflowmeter.
2. Place the Uroflowmeter on the charging puck to charge.
3. Ensure required items in Table 1 are obtained and set up.
4. Ensure mobile devices are connected to wireless internet and desktops are connected to internet.
5. Ensure users have proper access to the Glean Mobile App and Web App.
6. After steps 1 – 4 are complete, the GUS is ready to begin the Urodynamics procedure(s).

Bright Uro Equipment	Customer Provided Equipment
<ul style="list-style-type: none"> ✓ Uroflowmeter ✓ Charging Puck ✓ Power cables ✓ Sensor and Insertion Tool ✓ Uroflowmeter Quick Start Guide ✓ Owner's manual 	<ul style="list-style-type: none"> ✓ Laptop or computer for use with the Glean Web App ✓ Mobile device (tablet/phone) for use with the Glean Mobile App ✓ Urine collection cup ✓ Commode chair and funnel ✓ Materials required for aseptic insertion technique (such as a Foley Catheter Insertion Kit). ✓ Water-based lubricant and/or lidocaine gel. ✓ Biohazard bags

Table 1. Required Equipment

1.2.1 How Supplied

1.2.1.1 Sterility

The Sensor and Insertion Tool are provided STERILE (ethylene oxide [EO] sterilization). The sterile packaging should be inspected for visible damage prior to use. Do not use if damage is suspected. Do not reuse or attempt to re-sterilize.

1.2.1.2 Contents

GUS may be provided as two separate packages. One package contains the Sensor and Insertion Tool and a separate package contains the Uroflowmeter.

The Uroflowmeter package contains the following components:

- Uroflowmeter (ME Equipment, no applied parts) - IP54 rated
- Uroflowmeter Wireless Charger - IP54 rated
- Uroflowmeter Charger AC/DC Adapter (Manufacturer // PN: HDP Power// HDP12-MD-WUSB-4) (no ingress protection. Keep away from wet areas)
 - Input ratings: 90~264VAC, 47~63Hz, 12W max
 - Charger rated voltage, power: 5V, 5W

(The combination of the Uroflowmeter Charger and AC/DC Adapter make up the ME System)

The Sensor and Insertion Tool package contains the following components:

- Sensor (Type BF Applied Part)
- Insertion Tool (Sheath and Advancer)

The Glean Mobile App can be downloaded directly from the Google Play™ store for Android products and Apple App Store™ for iOS products. The Glean Web App can be accessed at gleanuds.com. The mobile apps can be downloaded directly from the Google Play™ store for Android products and Apple App Store™ for iOS products.

1.2.1.3 Additional Required Items

- Laptop or computer for use with the Glean Web App
- Mobile device (tablet/phone) for use with the Glean Mobile App
- Urine collection cup
- Commode chair and funnel
- Materials required for aseptic insertion technique (such as a Foley Catheter Insertion Kit)
- Water-based lubricant and/or lidocaine gel
- Biohazard bags

1.2.2 Device Inspection

Inspect each device and packaging to verify that no damage or defects exist. If the device is expired, damaged or if the sterile barrier has been compromised (e.g., hole in device packaging), do not use the device.

1.2.3 Device Storage

- Store the kits at room temperature.
- Avoid direct sunlight.

1.3 LEARNING ABOUT THE GLEAN URODYNAMICS SYSTEM (GUS)

To learn about the features and how-to of the Glean Urodynamics System, read the following documents or sections of this manual:

- GUS Uroflowmeter Quick Start Guide (provided with the equipment shipment)
- GUS Training Videos (may be accessed at gleanuds.com/training)
- Software Features and Functions (see 7 SOFTWARE FEATURES AND FUNCTIONS on page 16)
- How to Run Tests – CMG/PF Test, Uroflow Test, and Data Analysis (see 8 HOW TO RUN TESTS – CMG/PF TEST, UROFLOW TEST, AND DATA ANALYSIS on page 23)

1.4 INTENDED USE / INDICATIONS FOR USE

The Glean Urodynamics System (GUS) is a Urodynamic Analyzer System that is intended to quantify the pressure characteristics of the lower urinary tract. The system can be used to perform standard Urodynamic tests such as Uroflow, Cystometrogram (CMG), Urethral Pressure Profile (UPP), and Micturition Studies.

The major application of Urodynamics is the diagnosis of uncontrolled loss of urine (incontinence), abnormal urinary retention, or neurological cases of micturition disorder. The device is intended to be used as medical diagnostic equipment.

1.5 CONTRAINDICATIONS

- Use of GUS is contraindicated for any patient who is not a candidate for Urodynamic testing.
- The Sensor should not be used on patients who suffer from symptomatic urinary tract infections. Prior to testing, a urinalysis and urine culture should be considered to rule out the presence of infection.
- The Sensor should not be used on patients who suffer from a major stricture in the urethra.
- Single-use, disposable Sensors and Insertion Tools provided by Bright Uro are “sterile,” unless stated otherwise on the packaging label and instructions.

1.6 TARGET USERS

- Only technicians and clinicians trained in Urodynamics should operate this device. The operator must read the Owner’s Manual entirely and refer to any additional training materials before using the device.
- To reduce the potential for discomfort during the procedure and/or transient discomfort, dysuria and hematuria, technicians and physicians should explain any additional risks of the procedure to the patient.

- To reduce the risk of serious patient injury, it is vital that clinicians performing Urodynamics studies on patients with a Spinal Cord Injury be prepared to recognize and treat Autonomic Dysreflexia. Clinicians must monitor patients with Autonomic Dysreflexia for at least 2 hours after resolution of the episode.
- To reduce the risk of cross-contamination or urinary tract infection, clinicians should be knowledgeable and qualified in applying the appropriate aseptic technique during the intended use of the device. The use of prophylactic antibiotics is at the discretion of the clinician and the policies of the clinic/institution.
- To reduce the risk of serious patient injury, it is vital that technicians and clinicians performing Urodynamics studies be prepared to recognize and treat symptoms associated with vasovagal syncope (fainting) during Urodynamics procedures.

1.7 WARNINGS AND PRECAUTIONS

PRECAUTIONS

- ⚠ Bright Uro equipment and accessories are licensed by Governments, approved by Safety Agencies, and warranted to work only with each other.

CAUTION: UNITED STATES FEDERAL LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A LICENSED PHYSICIAN.

⚠ SYSTEM WARNINGS

- ⚠ DO NOT USE GUS in the presence of a magnetic resonance imaging (MRI) system as it may contain ferromagnetic objects that pose a risk to the patient in the presence of a magnetic core. The strong magnetic field produced by the MRI may cause disruption of the system.
- ⚠ DO NOT ATTEMPT TO OPEN OR REPAIR GUS components by yourself or by an unauthorized party. ONLY Bright Uro trained technicians may service GUS components.
- ⚠ Batteries are not operator removable; do not attempt to remove batteries from the GUS system components. All servicing of the GUS system, components or attachments are to be completed by Bright Uro.
- ⚠ Exposure to electrostatic discharge (ESD) may cause GUS to FAIL.
- ⚠ Bright Uro is not responsible for loss of patient files or test data.
- ⚠ Re-use, reprocessing or re-sterilization of disposables can lead to device failure and create a risk of cross-infection and/or cross transmission of infectious disease(s) from one patient to another. The Insertion Tool and Sensor are provided as single use, disposable devices and are intended to be discarded after use.
- ⚠ DO NOT immerse the GUS Uroflowmeter or other reusable system components in water or any other liquids. Only use approved cleaning agents to clean the Uroflowmeter as outlined in this Owner's Manual.
- ⚠ WARNING: Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- ⚠ WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Sensor or Uroflowmeter, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

ⓘ SYSTEM IMPORTANT INFORMATION

1. Use GUS with Bright Uro equipment and accessories only. Do not reuse disposable devices. After use, dispose in accordance with local regulations. Do not use if device packaging has been opened, or damaged, or if it presents any fault due to improper transport, storage, or handling that could in any way hamper its use.

2. Device intended for use in a clinical environment with controlled electromagnetic compatibility (EMC) standards to limit potential interference. GUS may be adversely affected by Bluetooth®, cellular or EMC interference. Minimize interference from other Bluetooth devices by setting up all components of system in proximity to each other. Placement of GUS Sensor on a patient's upper torso should be avoided to minimize any possibility of electromagnetic interference with active implantable devices such as ICD's and pacemakers.
3. The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
4. Frequency: The system radios operate in the 2.4 GHz ISM (Industrial, Scientific, Medical) band. This band specifically ranges from 2.400 to 2.4835 GHz. Bandwidth: The system radios use a frequency-hopping spread spectrum, where it occupies a channel with a bandwidth of 2MHz. Modulation: BLE primarily uses Gaussian Frequency Shift Keying (GFSK) modulation. Effective Radiated Power: +0dBm
5. The appliance coupler or AC/DC power adapter is not used as an isolation means. The Uroflowmeter isolation is reinforced through plastic and is decoupled to the charger. The 2 means of protection come from the plastics around the charger and the uroflowmeter. Also, the uroflowmeter does not have an applied part and does not contact the patient during normal use. The operator touches the device during use and charge. Therefore, the protection is means of operator, not patient.



SYSTEM SAFETY COMPLIANCE

1. To prevent unexpected exposure to radiation, the device has been tested against EN 60601-1-2 EMC standards.
2. To prevent exposure to potential electric shock, the device meets and exceeds the insulation breakdown specifications for IEC 60601-1:2005 & A2:2020; EN 60601-1:2006 & A2:2021; ANSI/AAMI ES60601-1:2005 & A2:2021; CAN/CSA C22.2 No. 60601-1:2014 & A2:2022.
3. To prevent radiofrequency electromagnetic interferences, the device meets and exceeds the specifications for IEC 60601-1-2:2014 & A1:2020; EN 60601-1-2:2015 & A2:2021
4. Warning symbols on all labels comply with ISO 7000, EN ISO 15223-1, and ISO 20417.
5. MRI Safety: As indicated on the device label, the Sensor is **MRI unsafe**. The Sensor should be removed before imaging or treatment.



NOTE: Local laws take priority over the above-mentioned requirements and warnings; if in doubt, consult your local Bright Uro representative or the technical service department.

2 ADVERSE EVENTS/RESIDUAL RISKS

2.1 POTENTIAL ADVERSE EVENTS/RESIDUAL RISKS

Possible complications associated with the use of GUS are similar to those associated with other methods of Urodynamics and include, but may not be limited to:

- Autonomic Dysreflexia (for individuals with spinal cord injury)
- Bladder or Urethral Spasms
- Change in Urinary Frequency
- Discomfort
- Injury to the Lower Urinary Tract (LUT) and/or Genitals
- Urinary Urgency
- Urinary Incontinence

- Urinary Retention
- Urinary Tract Infection
- Urinary Tract Inflammation or Irritation

2.2 DEVICE-RELATED ADVERSE EVENTS REPORTING

Any device-related adverse event or other incident regarding GUS should be immediately reported to Bright Uro. To report an event or incident, email: support@brighturo.com.

3 PATIENT COUNSELING INFORMATION

The physician should review the risks and benefits with the patient. Patients with a history of urethral strictures or frequent urinary tract infections should be monitored closely.

4 GLEAN URODYNAMICS SYSTEM—CARE AND MAINTENANCE

4.1 GENERAL CARE, CLEANING, AND PREVENTIVE MAINTENANCE

 **IMPORTANT:**

- Always wear protective gloves when cleaning the equipment to prevent biological contamination.
- GUS Sensors and Insertion Tools are intended for SINGLE PATIENT USE only. Do NOT reuse disposables.
- The Uroflowmeter component of GUS is non-immersible. The Uroflowmeter should be wiped down with a clean cloth dampened with a cleaning solution such as soap and water or as per hospital cleaning instructions.
- Do not sterilize the GUS components.
- Performing regular maintenance will reduce the need for costly repairs.
- Pay close attention to the LED lights on each device. If they indicate a broken connection and/or low battery, make sure that the connection is reestablished and fully charge the battery. Refer to 6.1 LED LIGHTS on page 16.

4.1.1 Caring for the GUS Uroflowmeter

The instructions that follow specify how to clean the GUS Uroflowmeter of possible urine contamination. Always wear protective gloves when cleaning the equipment to prevent biological contamination.

 **IMPORTANT:** Do not soak the GUS Uroflowmeter in water! Do not immerse in water or in any other liquids!

- The GUS Uroflowmeter has an immersion protection rating of IP54 for ingress of water. This means that the enclosure of the Uroflowmeter can handle splashes of water and liquid from any direction, but it is not protected against total immersion into water or any other liquids. This complies with IEC 60529 standards.
- The Glean Uroflowmeter should be cleaned using the following steps:
 - Separate the Uroflowmeter from the AC Adapter and Charging Cable.
 - Power off the Uroflowmeter by pressing and holding the Button LED on the front of the device.
 - Wipe the Uroflowmeter. Potential cleaning solutions in clinic include:
 - A cloth dampened with Isopropyl Alcohol
 - Soap
 - Disinfectant Detergent
- The Uroflowmeter may be used after it is wiped down.
- **WARNING: DO NOT SUBMERSE THE UROFLOWMETER IN ANY FLUID. DOING SO MAY DAMAGE OR DESTROY THE DEVICE.**
- Store the Glean Uroflowmeter in a cool and dry area at room temperature.
- Calibration: Return the Uroflowmeter annually to Bright Uro for recalibration. Contact Bright Uro to schedule this service as required.

4.2 BATTERY—CHARGING AND PREVENTATIVE MAINTENANCE

4.2.1 Charging the Battery

The GUS Uroflowmeter contains rechargeable batteries. The battery status is indicated by the Button LED on the front of the Uroflowmeter. For information on battery status, see 6 EQUIPMENT STATUS INDICATORS AND BUTTONS on page 16.

To charge the GUS Uroflowmeter:

- Plug the power cable of the Uroflowmeter Charging Puck into an electrical outlet.
- Place the GUS Uroflowmeter securely onto the Charging Puck.
- Confirm that the Uroflowmeter button LED shows that the device is being charged.



NOTE:

- When the Button LED on the device is BLUE, it is charging.
- The Uroflowmeter is the only rechargeable device in the Glean Urodynamics System.
- The Uroflowmeter can be used while charging.

4.4 TREATING AND DISPOSING OF PRODUCT AFTER USE

- After use, discard the contaminated, plastic, single-use disposables and any packaging according to your institution's standard operating procedures on medical-waste handling.
- For end-of-life product, waste electrical and electronic equipment should be collected separately and returned to the designated local recycling service.
- For end of battery life, disposal must be handled according to local regulations.
- Packaging waste should be collected separately for available national packaging collection and recycling services.

4.5 ENVIRONMENTAL CONSIDERATION OF WASTE DISPOSAL

Because the GUS is designed to perform Uroflow studies using the Uroflowmeter, it is important to dispose of waste (such as urine) properly to prevent environmental pollution. The waste should be disposed of in such a way that will not pollute the freshwater supply system—especially the drinking water system. In areas that have sewage systems with water treatment procedures, the most convenient method of disposal is to use these sewage systems.

4.6 PREVENTIVE MAINTENANCE—CHECKING CALIBRATION

Return the Uroflowmeter annually to Bright Uro for recalibration. Contact Bright Uro to schedule this service as required.

5 SET UP THE GLEAN URODYNAMICS SYSTEM (GUS)

Upon receiving the Glean Urodynamics System, inspect the equipment for any visible signs of damage or mishandling. If damage has been found, notify the carrier immediately. Bright Uro recommends saving carrying cases and cartons to provide a convenient and safe way to return the equipment should service be required.

The Uroflowmeter and its Charging Puck are the only reusable components of GUS. The following are required to set up the GUS Uroflowmeter for Urodynamics:

- Charge the Uroflowmeter by placing it on its Charging Puck.
- Place the Uroflowmeter on the floor as needed for tests.
- Press and hold the Button LED to power on the Uroflowmeter.
Place an unused urine collection cup on the Uroflowmeter. Ensure the urine collection cup is placed inside the silicon ring as indicated in Figure 5.
- To run a Uroflow Test refer to 8.1.9 Run a Uroflow Test on page 36



Figure 5. Urine Collection Cup Placement

The Sensor and Insertion Tool are intended for SINGLE PATIENT USE only. Do NOT reuse disposables. To set up the GUS Sensor for Urodynamics and to run a CMG/PF Test, refer to 8.1 CMG/PF TEST on page 23.

6 EQUIPMENT STATUS INDICATORS AND BUTTONS

6.1 LED LIGHTS

The LED light for the GUS Sensor is located in the middle of the unit. The LED lights for the GUS Uroflowmeter are located on the front of the device. Table 2 provides a description of LED light locations.




Table 2. Device LED Light Locations

6.1.1 Sensor LED Patterns

The Sensor LED will flash after the power button is pressed for 3 seconds while it is pairing via Bluetooth. Once paired, the Sensor will stop flashing. The Sensor will continue flashing for 150 seconds if it is not paired with a Bluetooth device.

6.1.2 Uroflowmeter LED Patterns

Device State	Ring LED	Button LED
OFF	OFF	OFF
ON	Pulsing	Pulsing
ON – Bluetooth Connected	Pulsing	Solid
ON – Collecting Data	Solid	Solid

Table 3. Uroflowmeter Device State LED Patterns

⚠ WARNING: If the Button LED displays a Triple RED Flash, then conditions have not been met to enter acquisition mode (e.g. battery state critical, Bluetooth not connected, session available, Power On Self-Test failed).

⚠ WARNING: If the Button LED displays a Single RED Flash, then the Power on Self-Test failed. Power the device off and power on again. If the problem persists, contact Bright Uro.

Battery State	Button LED
CHARGING	BLUE
NORMAL	WHITE
LOW	YELLOW
CRITICAL	RED

Table 4. Uroflowmeter Battery State LED Patterns

✔ NOTE: If the battery state is LOW or CRITICAL, place the Uroflowmeter on the charging puck.

7 SOFTWARE FEATURES AND FUNCTIONS

This section provides instructions on how to navigate the GUS Software – Glean Web App and Mobile Apps. For more information on the software, refer to 1.1.4 Software – Glean Mobile Apps (Clinician and Patient) and Glean Web App on page 9.

7.1 ACCESSING THE GLEAN APPS

Users can access the Glean Web App at gleanuds.com. The Glean Mobile App can be downloaded directly from the Google Play™ store for Android products and Apple App Store™ for iOS products.

7.1.1 Access the Glean Web App

Navigate to the Glean Web App at gleanuds.com

7.1.2 Install the Glean Mobile Apps

1. Login to mobile device.
2. Navigate to the Android/iOS App Store.
3. Search for the “Glean UDS” mobile app.
4. Download the Glean Mobile App.
5. Open the Glean Mobile App and allow the app to use Bluetooth.




Figure 6. Glean Mobile App Bluetooth Connection

6. Create an account or login using account information.

7.2 LOGIN TO THE GLEAN MOBILE OR WEB APPS

Users can log into the Glean Mobile or Web Apps with account information.

1. Open the Glean Mobile or Web Apps.
2. Enter user email.
3. Enter user password.
4. Click “Sign In” to login.




Figure 7. Glean Mobile App and Web App Login Page

 **NOTE:** If a device is provided to patients by the clinic, ensure the patient is logged out of account upon returning the device. Accounts will be automatically logged out after 30 minutes of inactivity.

7.3 CREATE USER ACCOUNTS

Clinic Admins (Referred to as Tenant in the Glean Web App) will create accounts for designated personnel to support GUS procedures.

1. Login to the Glean Web App (Admin).
2. Navigate to the desired account page (Admin/CP/Patient) using the three bars icon on the top left.




Figure 8. Navigate to Account Page

3. Select the three dots icon at the top right and click "CREATE."



Figure 9. Create an Account

4. Enter the required information and click "SUBMIT."




Figure 10. Enter Account User Information

7.4 DELETE USER ACCOUNTS

Clinic Admins may use the web portal to delete user accounts as required.

1. Login to the Glean Web App (Admin).
2. Navigate to the desired account page (Admin/HCP/Patient).




Figure 11. Navigate to Account Page

3. Locate the user account and click “VIEW.”




Figure 12. View User Account

4. Select the three dots icon at the top right and click “DELETE.”



Figure 13. Delete User Account

5. Click “YES” to delete user account.

Are you sure?

YES

NO

Figure 14. Confirm Deletion of User Account

7.5 RESET PASSWORD

Admin or users may reset a password for a Glean account.

7.5.1 Admin

1. Login to the Glean Web App (Admin).
2. Navigate to the desired account page (Admin/HCP/Patient).




Figure 15. Navigate to Account Page

3. Locate the user account and click "VIEW."




Figure 16. View User Account

4. Select the three dots icon at the top right and click "RESET PASSWORD."




Figure 17. Reset Password

- Observe the success message at bottom of screen saying, “Email successfully sent.”



Figure 18. Reset Password Email Successfully Sent

7.5.2 User

- Open the Glean Mobile or Web App.
- Select “Forgot Password?”




Figure 19. Forgot Password

- Enter user email and click “Send” or “Submit.”




Figure 20. Send Email to Reset Password

- Click on verification link in email.
- Enter a new password.
- Login to the Glean Mobile or Web App using the new password.

7.6 ADD A PATIENT PROFILE IN THE GLEAN ADMIN PORTAL

Use the Glean Admin portal to ensure the correct patient information is uploaded for future Urodynamics evaluation.

- Obtain patient information needed for entry into the Glean Web App (Admin).
- Login to the Glean Web App (Admin).

3. Select the three bars icon and click on “Patients.”




Figure 21. Patient Account Page

4. Select the three dots icon and click “CREATE.”



Figure 22. Create Patient Account

5. Enter the patient information and click “SUBMIT.”




Figure 23. Enter Patient Information

6. Observe the green pop-up window at the bottom of the screen to confirm patient entry.

 **NOTE:** If a patient’s email is already in the system but not associated with a clinic, a pop up will ask “Do you want to invite this patient to your clinic?” Select “OK”, and an email will be sent to the patient to ask if they would like to associate with the clinic. Have the patient accept the invitation to be associated with your clinic.

8 HOW TO RUN TESTS – CMG/PF TEST, UROFLOW TEST, AND DATA ANALYSIS

This section provides instructions on how to run tests with the GUS. For instructional videos on how to run tests with GUS, visit gleanuds.com/training. For more information on equipment or accessories setup, refer to 1 INTRODUCTION on page 8.

8.1 CMG/PF TEST

The purpose of running a CMG/PF test is to determine whether the bladder and its surrounding tissues are functioning correctly. The CMG test involves allowing the bladder to fill in a natural, antegrade manner and determining the vesical pressure, Pves, via the GUS Sensor.



NOTE: Make sure the batteries on the devices are fully charged before the start of the test.

8.1.1 Prepare the patient for aseptic insertion.

8.1.2 Instill lubricant in the urethra.

Instill lubricant with or without lidocaine in the urethra if needed based on clinical judgement.

8.1.3 Prepare the Sensor for data collection.

1. Inspect labeling to select the proper sensor for the patient's gender (male or female).
2. Open the outer box and remove pouch. Do not open the pouch or remove the Sensor from the pouch at this step.
3. Login to the Glean Mobile App (Clinician).
4. Select "Start a New Study."




Figure 24. Start a New Study

5. Select the correct patient.




Figure 25. Select a Patient

6. Confirm the patient's information.




Figure 26. Confirm Patient Information

7. Confirm the patient's information and select "Start Study."




Figure 27. Start Study

8. Select “Add Events or Devices.”




Figure 28. Add Events or Devices

9. Select “Add Sensor.”




Figure 29. Add Sensor

10. Press and hold the Sensor button for >3 seconds while inside the packaging to power on the sensor and click “Next” when complete. When the sensor is powered on you will see the LED flashing.




Figure 30. Power On the Sensor

11. Scan the QR code on the pouch and allow the Glean Mobile App to access your camera.




Figure 31. Scan the QR Code on the Pouch

12. Enter the PIN from the label near the QR Code. Once bonded, click “Setup Sensor.”





Figure 32. Enter PIN and Setup Sensor

13. Follow the instructions to calibrate the Sensor by leaving the Sensor on a flat surface. Click “Done” when calibration steps are complete.




Figure 33. Calibrate the Sensor

 **NOTE:** A pop up message will tell you if the Sensor is not functional and to discard. If this occurs, repeat the steps in 8.1.3 with a new Sensor.

8.1.4 Load the Sensor in the Sheath.

1. Peel the pouch open using the chevron end.
2. Gently drop the tray onto a flat surface.
3. Prepare for aseptic insertion (wash hands, put on sterile gloves, gather any additional items needed for a sterile field, etc.).
4. Open the plastic tray lid.
5. Apply lubricant to the entire length of the Sensor with specific focus on the space between the Sensor endcap and the Sheath.
6. Grab the Removal String and place in line with sheath so that it is hanging over the edge of the plastic tray.
7. Gently close the tray lid ensuring that each snap is properly seated, especially around the Sensor.
8. Using the left hand to hold the tray still, gently pull the Removal String until the sensor is fully loaded into the Sheath.
9. Open the plastic tray lid and lift the Sheath out of the plastic tray.
10. Inspect the Sensor tip to ensure proper alignment and seating of the Sensor and Sheath locking feature.
11. If necessary, apply gentle pressure to rotate and seat the Sensor into the Sheath locking feature.

8.1.5 Deploy the Sensor in the bladder.

8.1.5.1 Male Patient

1. Apply gentle traction to the tip of the penis to straighten the urethra.
2. Grasping the body of the Sheath, insert the tip of the Sensor into the urethra.
3. While maintaining traction on the penis, gently advance the Sheath ensuring not to push past any significant resistance.
4. Continue advancing the Sheath until the handle of the Sheath is near the tip of the penis.
5. Gently withdraw the Sheath approximately 2-4 centimeters.

6. Use one hand to grasp the handle of the Sheath and the other hand to pick up the Advancer.
7. Insert the Advancer and push to deploy the Sensor.
8. Continue pushing the Advancer until the handle meets the Sheath.
9. Gently withdraw the Advancer and confirm placement with visual observation of urine flow.
10. If urine does not flow, maintain the positioning of the Sheath and wait at least 20 seconds to observe urine flow.
11. If urine still does not flow, remove the Sheath then remove the sensor and reattempt insertion with a new Sensor once the patient's bladder has filled.
12. Gently remove the Sheath ensuring not to pull the Removal String.
13. Secure the Removal String to the patient's body using tape (or similar materials such as Tegaderm).
14. Click "Deployment Complete" or "Deployment Failure" on the Glean Mobile App (Clinician) when Sensor deployment is completed/failed (Figure 34).

8.1.5.2 Female Patient

1. Separate the labia to expose the urethra.
2. Grasping the body of the Sheath, insert the tip of the Sensor into the urethra.
3. Continue advancing the Sheath until approximately ½ of the Sheath is inside the patient's body.
4. Use one hand to grasp the handle of the Sheath and the other hand to pick up the Advancer.
5. Insert the Advancer and push to deploy the Sensor.
6. Continue pushing the Advancer until the handle meets the Sheath.
7. Gently withdraw the Advancer and confirm placement with visual observation of urine flow.
8. If urine does not flow, maintain the positioning of the Sheath and wait at least 20 seconds to observe urine flow.
9. If urine still does not flow, remove the Sheath then remove the sensor and reattempt insertion with a new Sensor once the patient's bladder has filled.
10. Gently remove the Sheath ensuring not to pull the Removal String.
11. Secure the Removal String to the patient's body using tape (or similar materials such as Tegaderm).
12. Click "Deployment Complete" or "Deployment Failure" on the Glean Mobile App (Clinician) when Sensor deployment is completed/failed (Figure 34).



NOTE: If at any time you feel resistance do NOT force the Insertion Tool. You may need to apply more lubrication before continuing insertion.




Figure 34. Sensor Deployment Complete or Failure

8.1.6 Log events using the Glean Mobile App (Clinician).

1. Login to the Glean Mobile App (if not already logged in).

- 2. If desired, perform any series of guided maneuvers based on patient history, symptom presentation, and goals for urodynamic evaluation.
- 3. Use the Glean Mobile App (Clinician) to select the correct patient for the ongoing study and select “Add Events or Devices.”




Figure 35. Add Events or Devices

- 4. Select “Add Events.”




Figure 36. Add Events (Clinician)

- 5. Select the desired event (Cough, Valsalva, PVR, 1st Sensation, 1st Desire, Strong Desire, Voiding Event, Leak, Fluid Intake, Urge, Others) and guide the patient to perform maneuver if necessary.





Figure 37. Select Appropriate Event (Clinician)

- 6. Confirm the event details and enter required information. Select “Save.” Examples of 1st Sensation and Leak are shown below in Figure 37.







Figure 38. Confirm Event Details (Clinician)

7. To log a Cough or Valsalva, ask the patient to perform maneuver after you click “Start Timer.” Click “Stop Timer” when patient is finished performing maneuver. An example of Cough is shown below in Figure 38.






Figure 39. Event Timer

8. For all other maneuvers, attempt to log the event and press the event button at the same time as when the patient is performing the maneuver.
9. Repeat steps 7-8 for each maneuver that is necessary given the patient’s history and symptom presentation.
10. Edit the date and time of each event by selecting an event from the Ongoing Study page and then selecting the pen icon. An example for editing the date and time is shown in Figure 40. After selecting the pen icon, you can modify the time or other data fields and save the changes. To do this, it will require a comment to justify the change to the event.
11. To delete an event, hold your finger on the event for at least one second and select “yes” to delete the event.

8.1.7 Instruct the patient to log symptoms during ambulatory monitoring.

1. If desired, have the patient log symptoms during the ambulatory monitoring period.
2. Ensure the patient has downloaded the Glean Mobile App to their smartphone or provide the patient with a device that has the Glean Mobile App installed. If preferred, the patient may use a pen and paper to log symptoms. A Symptom Log template is available at gleanuds.com/diary.
3. Ensure the patient has logged into the Glean Mobile App with the proper account information.
4. Educate the patient on how to log events (refer to Chapter 11.1.8 on page 27 for details on how to log events using the Glean Mobile App (Patient)).




Figure 40. Edit Event (Clinician)

5. Instruct the patient to return to the exam room when they have a strong desire to void.

8.1.8 Train the Patient to use the Glean Mobile App (Patient).

1. Login to the Glean Mobile App (Patient).
2. Select "Add Event."




Figure 41. Add Event (Patient)

3. Select event type (leak, urge, fluid intake, voiding event) and enter data.




Figure 42. Select Event Type (Patient)

4. Enter required information, Select “Enter” then “Save” to complete data upload.




Figure 43. Enter Event Details (Patient)

5. Edit a logged event by selecting an event from the home page and clicking “Edit” or by clicking on the pen icon when entering the data.




Figure 44. Edit Event (Patient)

6. Delete a logged event by clicking “Delete.”




Figure 45. Delete Event (Patient)

8.1.9 Run a Uroflow Test

A Uroflow Test is a measurement of the rate at which urine flows out of the body. It can be performed using the Glean Mobile App (Clinician). A Uroflow Test may be conducted as a standalone test, when a Sensor is not in use, or as part of a CMG/PF Test when a Sensor has been deployed in the patient.

 NOTE: Make sure the battery of the GUS Uroflowmeter is fully charged before starting the test.

8.1.9.1 Prepare the Uroflowmeter for data collection.

1. Gather the supplies needed for a Uroflow Test (urine collection cup, commode chair, funnel, etc.).
2. Carefully place the GUS Uroflowmeter on the floor.
Gently position a urine collection cup on top of the Uroflowmeter. Ensure the urine collection cup is placed as indicated in
3. Figure 5.
4. Place the funnel on the plastic frame of the commode chair and position both over the Uroflowmeter and receptacle. Ensure that the urine collection cup and the funnel are aligned, but not touching.
5. Login to the Glean Mobile App (Clinician).
6. If adding Uroflowmetry to an existing study, select the ongoing study from the active studies list and go to Step 9.

7. If required, select “Start New Study.”




Figure 46. Start a New Study

8. Select the desired patient profile, confirm the patient information, and start the study if necessary.




Figure 47. Select a Patient, Confirm Patient Information, and Start Study

9. Select “Add Events or Devices.”




Figure 48. Add Events or Devices

10. Select “Add Uroflowmeter.”




Figure 49. Add Uroflowmeter

11. If required, power on the Uroflowmeter by pressing and holding the button for at least 3 seconds and click “OK” when complete.




Figure 50. Power on the Uroflowmeter

12. Scan the QR code on the Uroflowmeter or Quick Start Guide.



Figure 51. Scan QR Code on Uroflowmeter

13. If required, enter the PIN near the QR code.




Figure 52. Enter Uroflowmeter PIN

14. Click "OK" to associate the Uroflowmeter with the patient profile.




Figure 53. Associate Uroflowmeter with Patient

15. Click “Connect” then ‘Start Uroflowmetry’ to start collecting Uroflowmeter data.




Figure 54. Start Uroflowmetry

 **NOTE:** If you lose Bluetooth connection during the void then you will need to repeat steps 5 – 14 after the void is complete to download the data to the patient profile, which may include stopping the study.

8.1.9.2 Instruct the patient to void.

1. Instruct the patient not to touch or kick the urine collection cup before, during or after voiding.
2. Tell the patient to void into the collection cup until they feel their bladder is empty and to notify clinic staff once complete.

 **CAUTION:** DO NOT TOUCH the urine collection cup during voiding.

8.1.9.3 Download the data from the Uroflowmeter.

1. Click “Stop Uroflowmetry” to stop data acquisition from the Uroflowmeter.




Figure 55. Stop Uroflowmetry

2. Click "Download" to download the Uroflowmetry data.




Figure 56. Download Uroflowmetry Data

3. View Uroflowmetry data and select “Confirm” to upload to the patient record.




Figure 57. Confirm Uroflowmetry Data

8.1.10 Measure PVR using preferred method.

1. If desired, measure PVR using preferred method.
2. Log bladder volume as an event titled “PVR” in the Glean Mobile App (Clinician).

8.1.11 Prepare the patient for Sensor removal.

1. Have the patient remove clothing.
2. Have the patient assume a comfortable position for Sensor removal.

8.1.12 Remove the Sensor from the bladder.

1. Gently remove any material used to secure the Sensor Removal String.
2. Gently pull the removal string until the Sensor is completely out of the body.
3. Place the Sensor in a biohazard bag and close the bag.

8.1.13 Urethral Pressure Profile (UPP)

If desired, Urethral Pressure Profile testing may be performed with Glean using a Manual Pull.

1. Ensure the patient has at least 50 mL of urine in the bladder.
2. Gently remove any material used to secure the Sensor Removal String.
3. Pull the Removal String very gently until you begin to feel resistance at the bladder neck.
4. Once you feel resistance from the Sensor at the bladder neck begin pulling very slowly. Continue to watch the Sensor as you remove it from the body and attempt to pull at a rate of approximately 1 mm per second.
5. Place the Sensor in a biohazard bag and close the bag.

8.1.14 Download data from the Glean sensor.

1. Select the correct patient profile from Glean Mobile App (Clinician) for the ongoing study.

2. Select "Add Events or Devices"




Figure 58. Add Events or Devices

3. Select "Download Sensor Data"




Figure 59. Select Download Sensor Data

4. Press and hold the sensor button for at least three seconds until the LED begins to flash.
5. Select "Connect."





Figure 60. Connect Sensor to Download Data

6. Select "Download" and keep the Sensor near the mobile device until the data download is complete.






Figure 61. Download Sensor Data

7. View Sensor data and select "Confirm" to upload data to the patient record.




Figure 62. Confirm Sensor Data

8. Dispose of the Sensor according to clinic guidelines.

8.2 ANALYZE THE DATA USING THE GLEAN WEB APP

Health Care Providers may analyze GUS data using the Glean Web App.

8.2.1 Login to the Glean Web App (Clinician).

8.2.2 Select the desired patient and study.

1. Current or pending studies will appear under the Studies Overview page in the default homepage view. Select the desired patient and study by clicking on the desired row

Status	MRN	First Name	Last Name	Study Start Date	Study Type	Last Activity	Time of Last Activity
Current	2670408-9370-4165-9267-2762...	Melissa	McAlpine	Mar 19, 2024 10:17:31	Cystometry	Volative	Mar 19, 2024 10:34:44
Pending	2670408-9370-4165-9267-2762...	Melissa	McAlpine	Mar 18, 2024 14:09:02	Cystometry	Volative	Mar 19, 2024 10:34:44
Pending	2222222	Bryan	Newspool	Mar 18, 2024 12:46:11	Cystometry	HearPressure	Mar 18, 2024 12:53:50
Current	12345555	teahyrz	vytest	Mar 14, 2024 02:58:29	Cystometry		
Pending	8770666-3529-4869-9c66-644278c05a8f	William	Bright	Mar 13, 2024 08:54:37	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	2222222	Bryan	Newspool	Mar 11, 2024 15:12:29	Cystometry	HearPressure	Mar 18, 2024 12:53:50
Pending	2670408-9370-4165-9267-2762...	Melissa	McAlpine	Mar 11, 2024 15:12:29	Cystometry	Volative	Mar 19, 2024 10:34:44
Pending	2222222	Bryan	Newspool	Mar 11, 2024 14:45:28	Cystometry	HearPressure	Mar 18, 2024 12:53:50
Pending	2670408-9370-4165-9267-2762...	Melissa	McAlpine	Mar 11, 2024 10:53:45	Cystometry	Volative	Mar 19, 2024 10:34:44
Pending	8770666-3529-4869-9c66-644278c05a8f	William	Bright	Mar 11, 2024 10:11:00	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22

Figure 63. Studies Overview

2. To view all patients, select the patient icon  on the left to view the list of patients.
3. Identify the correct patient using name or MRN and click "Studies."

MRN	First Name	Last Name	Status	Last Study Type
032942	Patient	82	Active	
test_patient_533	Robert	Law	Active	cystometry
517	Anthony	Shark	Active	cystometry
test_patient_343@example.com	test_patient_343	test_patient_343	Active	cystometry
test_patient_289@example.com	test_patient_289	test_patient_289	Active	cystometry
test_patient_254@example.com	test_patient_254	test_patient_254	Active	cystometry
test_patient_270@example.com	test_patient_270	test_patient_270	Active	
test_patient_888@example.com	test_patient_888	test_patient_888	Active	cystometry
test_patient_182@example.com	test_patient_182	test_patient_182	Active	
test_patient_182@example.com	test_patient_182a	test_patient_182a	Active	

Figure 64. Select Patient Studies

4. Confirm the date/time of the study and the type of study then click on the desired row.

Status	MRN	Study Start Date	Study Type	Last Activity	Time of Last Activity
Current	8770666-3529-4869-9c66-644278c05a8f	Mar 1, 2024 14:26:36	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Mar 1, 2024 09:30:35	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Feb 19, 2024 07:07:57	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Reviewed	8770666-3529-4869-9c66-644278c05a8f	Feb 9, 2024 09:05:42	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Jan 26, 2024 11:05:30	Uroflowmetry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Jan 26, 2024 11:01:49	Uroflowmetry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Jan 23, 2024 13:49:13	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Jan 23, 2024 13:32:17	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Jan 19, 2024 11:16:16	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22
Pending	8770666-3529-4869-9c66-644278c05a8f	Jan 19, 2024 09:05:01	Cystometry	VoidingEvent	Mar 5, 2024 12:01:22

Figure 65. View Patient Studies

5. The study data will open to be viewed.



Figure 66. View Study Data

8.2.3 Analyze Urodynamics data.

1. Visually review Urodynamics data.
2. Adjust or modify view as required for detailed analysis.

8.2.4 Adjust the x-axis for Urodynamics data.

1. Determine area of interest.
2. Select the + Zoom icon  at the top right of the viewing area.
3. Drag the box over desired area of interest to zoom in.




Figure 67. Drag Zoom Box

OR

1. Determine area of interest.
2. Adjust the left side of sliding view window to set the left limit for desired area of interest.

3. Adjust the right side of sliding view window to set the right limit for desired area of interest.



Figure 68. Adjust Sliding View Window

OR

1. Identify event of interest (e.g. cough, leak, void).




Figure 69. Select Event of Interest

2. Click on event of interest. The view will automatically zoom to center event with buffer on each side.




Figure 70. View Event of Interest

8.2.5 Adjust the y-axis for Urodynamics data.

1. Determine area of interest.
2. Adjust x-axis to zoom in to area of interest.
See Chapter 11.2.4 on page 41 for instructions on how to adjust the x-axis.
3. Software will automatically adjust to minimum and maximum pressure levels in area of interest.

8.2.6 Select multiple events of similar type.

1. Navigate to the Uroflowmetry or Voiding Diary tab.



Figure 71. Uroflowmetry and Voiding Diary Tabs

2. Select the desired events from the table in the bottom viewing area.

UF	Start time	Voided Volume (fl oz)	Qmax (ml/s)	T-Qmax (sec)	Flow Time (sec)	End time
<input checked="" type="checkbox"/> UF #avg	13:59:41	9.58	0.44	151.35	-90.04	14:03:35
<input checked="" type="checkbox"/> UF #0	13:59:41	9.58	0.44	151.35	-90.04	14:03:35

<input checked="" type="checkbox"/>	Date	Fluid Intake (fl oz)	Volume Voided (fl oz)	Voiding Events	Nocturia Events	Pad Use	Urgency Severity	Leakage Events
<input checked="" type="checkbox"/>	2024-01-19	8.00	0.00	0	0	0	0	1
<input checked="" type="checkbox"/>	2024-01-23	0.00	285.00	1	0	0	0	0

Figure 72. Select Events of Similar Type

3. Unselect events to be removed from the visual display and analysis.

UF	Start time	Voided Volume (fl oz)	Qmax (ml/s)	T-Qmax (sec)	Flow Time (sec)	End time
<input type="checkbox"/> UF #avg	13:59:41	9.58	0.44	151.35	-90.04	14:03:35
<input checked="" type="checkbox"/> UF #0	13:59:41	9.58	0.44	151.35	-90.04	14:03:35

<input type="checkbox"/>	Date	Fluid Intake (fl oz)	Volume Voided (fl oz)	Voiding Events	Nocturia Events	Pad Use	Urgency Severity	Leakage Events
<input type="checkbox"/>	2024-01-19	8.00	0.00	0	0	0	0	1
<input checked="" type="checkbox"/>	2024-01-23	0.00	285.00	1	0	0	0	0

Figure 73. Unselect Events of Similar Type

8.2.7 Draft notes for interpretation, assessment, and plan.

1. Type notes in text window for “Primary interpretation.”



Figure 74. Draft Primary Interpretation

2. Type notes in text window for “Assessment and plan.”




Figure 75. Draft Assessment and Plan

3. Select the copy icon to copy text for use outside of the Glean Web App.




Figure 76. Copy Notes

4. To save the notes, select the save icon or click outside of the notes text box to autosave.




Figure 77. Save Notes

8.2.8 Export the Urodynamics report.

1. Complete analysis and notes for the Urodynamics evaluation.
2. Confirm report accuracy.

3. Select "Complete Report."




Figure 78. Complete Uroynamics Report

4. Select "View Report."




Figure 79. View Uroynamics Report

5. Select "Print Report."




Figure 80. Print Uroynamics Report

6. Select desired destination for printing of report or “Save as PDF” and confirm file location.

8.2.9 Utilize filtering to adjust view of Urodynamics data.

1. Click on “Filters” in the top left.




Figure 81. View Filters

2. Determine desired events or desired filtering methods.




Figure 82. Select Filters

3. Select desired events/filters to display data in view.

4. Unselect event/filters to hide data from view.

8.2.10 View the Glean patient history.

1. Login to the Glean Web App (Clinician).

2. Select a patient record and study.




Figure 83. Select Patient Record and Study

3. Select "View Patient Glean History."



Figure 84. View Glean Patient History

8.2.11 View detailed Urodynamic study data.

1. Login to the Glean Web App (Clinician).
2. Select a patient record and study.




Figure 85. Select Patient Record and Study

3. Select “Study Data” to see detailed parameters of Urodynamics and Uroflowmetry exams.




Figure 86. View Study Data

8.2.12 Adjust Glean Web App settings.

Users may select or unselect settings to adjust the view of urodynamics data.

1. Login to the Glean Web App (Clinician).
2. Select the user icon  at the top right of the screen and click on “Settings.”




Figure 87. Glean Web App Settings

3. Adjust the user settings to accommodate the preference of the user and click “Save Settings.”




Figure 88. Adjust Glean Web App Settings

4. Click “Yes” to confirm your changes.




Figure 89. Confirm Settings Changes

9 CALIBRATION

Return the Uroflowmeter annually to Bright Uro for recalibration. Contact Bright Uro to schedule this service as required.

DRAFT

10 FAQ

How do I know when I need to charge the Uroflowmeter?

The GUS Uroflowmeter's LED lights indicate equipment status. If the Button LED is yellow or red, charge the Uroflowmeter. The Uroflowmeter is charging when the Button LED turns blue.

Can I use the Uroflowmeter while the battery is charging?

Yes, the GUS Uroflowmeter can be used while charging.

What should I do if the Removal String falls out of the Insertion Tool prior to loading the Sensor into the Insertion Tool?

Use graspers or forceps to pull the Removal String through the Sheath while maintaining a sterile field or use a new Sensor.

Can I use a dilator to assist in inserting the Sensor?

Yes, you may use a dilator to assist in inserting the Sensor if necessary, based on clinical judgement.

How do I know if the Sensor is properly deployed in the bladder?

The Sensor is properly deployed in the bladder if urine flows through the Sheath after Sensor deployment.

Is the Sensor transmitting Bluetooth data through the body?

The Sensor only connects via Bluetooth before insertion and after removal. The data is stored on the Sensor while it is indwelling in the body.

How do I dispose of the Sensor?

The Sensor is a single use device. Dispose of the Sensor according to local guidelines.

Are mobile devices and desktops provided by Bright Uro?

Bright Uro does not provide mobile devices or desktops. Clinics can utilize any available mobile devices and desktops.

11 TROUBLESHOOT

Symptom(s)	Possible Cause(s)	Check/Corrective Action(s)
SENSOR AND INSERTION TOOL		
Sensor LED is not flashing after pressing and holding power button for at least 3 seconds	Power button not pressed firmly	Firmly press and hold the power button for at least 3 seconds. Ensure proper finger placement on the power button.
	Sensor battery died	Use a new Sensor.
Unable to load the Sensor in the Insertion Tool	Removal String broken	Use a new Sensor.
	Sensor not aligned with locking feature	Twist the Insertion Tool while maintaining aseptic technique to align the locking feature with the Sensor.
	Lubricant not applied to Sensor	Apply lubricant over the entire Sensor and continue loading.
	Lid not closed over Sensor and Insertion Tool prior to loading	Close lid over Sensor and Insertion Tool and continue loading.
Unable to insert device into bladder	Sheath met with resistance	Gently advance the Sheath during insertion. Do NOT push past significant resistance. Reattempt or stop device insertion based on clinical judgement.
Unable to deploy Sensor in bladder	Urine did not flow through the Sheath	Maintain the positioning of the Sheath and wait at least 20 seconds to

		observe urine flow. If urine still does not flow, remove the Sheath then remove the sensor and reattempt insertion with a new Sensor once the patient's bladder has filled.
	Advancer met with resistance	Gently advance the Advancer during Sensor deployment. DO NOT push past significant resistance. Withdraw both the Sheath and Advancer together 2-3 cm then continue deployment. Reattempt or stop device insertion based on clinical judgment.
Unable to connect Sensor to Glean Mobile App (Clinician) to start CMG/PF Test.	Sensor LED not flashing	Press and hold the power button for at least 3 seconds until the LED begins to flash.
	Unable to scan QR code	Lay the Sensor packaging on a flat surface. Ensure the QR code is completely visible and rescan.
	Unable to connect via Bluetooth	Reattempt to connect. If still unable to connect, use a new Sensor.
Unable to connect Sensor to Glean Mobile App (Clinician) to download data.	Sensor LED not flashing	Press and hold the power button for at least 3 seconds until the LED begins to flash.
	Unable to connect via Bluetooth	Reattempt to connect. If still unable to connect contact Bright Uro.
UROFLOWMETER		
Uroflowmeter will not power on	Uroflowmeter not charged.	Place the Uroflowmeter on the Charging Puck and plug the power cable into an electrical outlet.
	Button LED not pressed firmly.	Press and hold the Button LED on the front of the device for at least 3 seconds.
Abnormal LED patterns	Conditions have not been met to enter acquisition mode	Ensure the device is fully charged. Power off and on the Uroflowmeter and reconnect.
	Power On Self-Test failed	Power off and on the Uroflowmeter again.
Unable to connect Uroflowmeter to Glean Mobile App (Clinician) to start Uroflow Test	Unable to scan QR code	Place the device on a flat surface. Ensure the QR code is completely visible and rescan.
	Unable to connect via Bluetooth	Forget device in Bluetooth setting of mobile device then reconnect with PIN.
	Uroflowmeter has previous data pending download.	Download data to correct patient profile and reattempt connection for new test.
	Uroflowmeter not powered on.	Press and hold the Button LED on the front of the device to power on.
	Uroflowmeter not charged.	Place the Uroflowmeter on the Charging Puck and plug the power cable into an electrical outlet.
Unable to connect Uroflowmeter to Glean Mobile App (Clinician) to download data	Uroflowmeter associated with different patient.	Ensure the Uroflowmeter being used is associated with the correct patient.

	Uroflow Study has not been stopped.	Stop study then download data.
	Uroflowmeter not powered on.	Press and hold the Button LED on the front of the device to power on.
	Uroflowmeter not charged.	Place the Uroflowmeter on the Charging Puck and plug the power cable into an electrical outlet.
SOFTWARE		
Unable to log in to Glean Mobile App or Web App	Incorrect password	Reset password and reattempt log in. If still unable to log in, contact clinic admin.
Unable to access the Glean Mobile or Web Apps	Unable to open Glean Mobile App	Power off and on mobile device or delete and redownload the app.
	Incorrect web address	Ensure navigation to correct web address via log in page located at gleanuds.com.

Table 5. Troubleshooting

Error Message	Location	Corrective Action(s)
WEB / MOBILE APP		
Could not get your token	Landing page	
MOBILE APP		
Connection failed	Sensor Setup page	Reattempt to connect. If still unable to connect contact Bright Uro.
There was an issue uploading sensor data to the cloud.	Data Download page	Wait for 2-3 minutes and reattempt to upload the sensor data to the cloud. If unsuccessful after 3 attempts, contact the Bright Uro Service team.
Error trying to read the QR Code	Scan Device QR Code	Make sure there is good lighting, the QR code symbol is completely visible, and the camera lens is clean. Reattempt to read the QR code
Unknown issue.	Sensor Setup page	Record the steps taken to get the error and contact the Bright Uro Service team.
Could not remove the Event	Edit Event page	Log out of the current session. Close the Glean UDS App and login to start a new session, then reattempt to remove the event. If the error persists, check the troubleshooting section or contact the Bright Uro Service team.
Error trying to save information	Add Event page	Log out of the current session. Close the Glean UDS App, and login to start a new session. Then reattempt to add the event. If the error persists, check the troubleshooting section or contact the Bright Uro Service team.
Error trying to delete your account	Patient Profile page	Contact the Bright Uro Service team for a possible solution.
Error trying to make the request.	Reset Password page	Log out from the current active session and close the Glean UDS App. Start a new browser and login to start a new session. Reattempt to reset the password. If the error persists, check

		the troubleshooting section or contact the Bright Uro Service team.
Error getting Clinics	Select Clinic page	
You must fill your Email.	Sign In/ Login page	Verify the correct email is entered.
You must provide your password.	Sign In/ Login page	Make sure the password field is filled out.
Limited or no internet connectivity.	Sign In/ Login page	Make sure you have an active internet connection and try to login again.
One event already has this name.	New Event Type (Others)	Use different event name that has not been used before and try again.
You cannot create an empty field event	New Event Type (Others)	Make sure all the required fields are filled out prior to creating the event.
You need to choose a Type	New Event Type (Others)	Make sure the event type is filled out.
Email address cannot exceed 64 characters after the "@"	New Event Type (Others)	Make sure to use a valid email address.
Email exceeds maximum length	New Event Type (Others)	Make sure to use a valid email address.
The selected hour cannot be in the future	New Event Type (Others)	Make sure to enter the valid time.

Table 6. Error Messages

If problems continue, contact the BRIGHT URO Service team at +1 (949) 216-0873 or by email at support@brighturo.com.

wireless

12 APPENDICES

APPENDIX A: TECHNICAL DATA

Model	Glean Urodynamics System
--------------	--------------------------

Classification EN 60601-1	- Applied part, Type BF - IP54 Rated
Mode	Continuous Operation
Sterilization	Ethylene Oxide (EO) Sterilization (Sensor and Insertion Tool only)
Operating Conditions	+10 °C (50 °F) to +40 °C (104 °F)(Sensor and Insertion Tool) +10 °C (50 °F) to +40 °C (104 °F) (Uroflowmeter) 15 % to 90 % Relative Humidity, non-condensing, but not requiring a water vapor partial pressure greater than 50 hPa (all devices)
Operating Atmospheric Pressure	700 - 1060 hPa
Transport and Storage Conditions	0 °C (32 °F) to +60 °C (140 °F) (Sensor and Insertion Tool, storage and transport) -20 °C (-4 °F) to +50 °C (122 °F) (Uroflowmeter, storage) -20 °C (-4 °F) to +70 °C (158 °F) (Uroflowmeter, transport) 5 % to 90 % Relative Humidity, non-condensing, but not requiring a water vapor partial pressure greater than 50 hPa, and 100 hPa to 1,060 hPa (Sensor and Insertion Tool) 5 % to 95 % Relative Humidity, non-condensing (Uroflowmeter, storage and transport) Rechargeable batter(ies) in Uroflowmeter to be removed prior to shipment or storage outside of the operational range. NOTE: Manufacturer considers that there is no hazard if device is used immediately after storage.
Weight	0.3 lbs (150 g)
Dimensions (H X W X D)	Sensor and Insertion Tool: 4.5" (114 mm) H x 3.5" (88 mm) W x 1" (25 mm) D Uroflowmeter: 3" (76.2 mm) H x 7.5" (190.5 mm) W x 7.5 (190.5 mm) D

Table 7. GUS Specifications

APPENDIX B: CLASSIFICATIONS

IEC 60601-1:	Class I Equipment Type BF Applied Parts
Mode of Operation:	Continuous; Equipment not suitable for use in the presence of a Flammable, Anaesthetic Mixture, with Air, or Oxygen, or Nitrous Oxide.
Degree of Protection:	The GUS Uroflowmeter enclosure is classified IP54 according to degree of protection against ingress of water and particulate matter as per the test requirements of IEC 60529. With this IP (International Protection) rating, it means that the Uroflowmeter enclosure: <ul style="list-style-type: none"> • Protects users using tools 1.0 mm or larger from accessing hazardous parts, and protects equipment from ingress of dust (signified by the rating code 5). • Protects equipment from the harmful effects of water splashing from any direction (signified by the rating code 4). NOTE: the IP rating will be visible on the Uroflowmeter label.

Table 8. GUS Classifications

APPENDIX C: SYMBOLS AND LABELING

	Consult Instructions for Use	STERILE EO	Sterilized using ethylene oxide
	Warning		Non-Sterile
	Do not use if package is damaged		Date of Manufacture
	Consult Instructions for Use		Manufacturer
	Not made with natural rubber latex	EC REP	Authorized Representative in the European Community
	Consult the instructions for use for important cautionary information such as warnings and precautions	Contents	Contents of box or container
	Strong Magnetic field	QTY	Quantity
	Keep dry	LOT	Lot number of product
	Keep out of direct sunlight	REF	Catalog number of product
	Single use only	SN	Serial number of product
	Humidity limitation		Use by Date (Expiration date of product)
	Upper Limit of Temperature		Atmospheric pressure limitation
	Temperature limits		Requires Disposal per Waste Electrical and Electronic Equipment Directive
	Type BF Applied Part		Peel open pouch at marked corner
R_x only	Prescription use only		Direct Current
	SGS Certification Mark to signify certification to applicable US safety standards		

1. ISO 15223-1 Medical Devices – Symbols to be used with medical device, labels, labelling and information to be supplied – Part 1: General Requirements.
 2. ISO 7010 Third Edition 2019-07 Graphical symbols – Safety colors and safety signs – Registered safety signs
 3. ISO 7000 Sixth edition 2019-07 Graphical Symbols for Use on Equipment – Registered Symbols.
 4. IEC 60417 – Graphic Symbols for Use on Equipment.
- NOTE: Sterility symbols are applicable to single-use devices only.

Table 9. Symbols Glossary

Labels can be found on the devices as follows:

Label Description	Label Placement
GUS Sensor and Insertion Tool	
GUS Uroflowmeter	

Table 10. Label Location

APPENDIX D: ELECTROMAGNETIC COMPATIBILITY (EMC)

This equipment has been tested and found to comply with the limits for:
 IEC 60601-1-2:2020(AMD+1), IEC 60601-2:40:2016 Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests

CISPR 11	Industrial, scientific, and medical (ISM) radio-frequency equipment – Electromagnetic disturbance characteristics
IEC 61000-3-2	Limits for harmonic current emissions (equipment input current = 16 A per phase)
IEC 61000-3-3	Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current = 16 A per phase and not subject to conditional connection
IEC 61000-4-2	Testing and measurement techniques – Electrostatic discharge immunity test

IEC 61000-4-3	Testing and measurement techniques – Radiated, radiofrequency, electromagnetic field immunity test. Ed 3.2.
IEC 61000-4-39	Testing and measurement techniques – FID Magnetic Proximity Fields
IEC 61000-4-4	Testing and measurement techniques – Electrical fast transient/burst immunity test
IEC 61000-4-5	Testing and measurement techniques – Surge immunity test
IEC 61000-4-6	Testing and Measurement Techniques – Immunity to Conducted Disturbances, Induced by Radio-Frequency Fields.
IEC 61000-4-8	Testing and measurement techniques – Power frequency magnetic field immunity test
IEC 61000-4-11	Testing and measurement techniques – Voltage dips, short interruptions, and voltage variations immunity tests
Clause 5	Identification, Marking and Documents

Table 11. Electromagnetic Compatibility

- These limits are designed to provide reasonable protection against harmful electromagnetic or other interference in most installations. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful electromagnetic or other interference, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate GUS Uroflowmeter unit.
 - Increase the separation between GUS Uroflowmeter unit and the affected equipment.
 - Connect the non-medical system equipment into an outlet on a circuit different from that to which the GUS Uroflowmeter unit is connected.
 - Consult experienced technical personnel for help.

 **WARNING:** Changes or modifications not expressly approved by Bright Uro could void the user’s authority to operate the equipment.

- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference received, including interference that may cause undesired operation.
- This device contains: FCC ID: PVH0946 IC: 5325A-0946

IEC 60601-1-2:2020 Table 1 Requirements

The GUS is intended for use in the electromagnetic environment specified below. The customer or the user of the GUS should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The GUS uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The GUS meets the conducted and radiated performance requirements for non-life supporting equipment and meets the harmonic emissions, voltage dips and variations and voltage fluctuation (flicker) requirements for non-life supporting equipment pursuant to CISPR 11, AI & A2, and IEC 61000-3-3.
Harmonic Emissions IEC 61000-3-2	Class A	

Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	<p>The GUS is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:</p> <p>Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the device or shielding the location.</p>
---	----------	--

Table 12. Table 1 Requirements—Electromagnetic Environment

IEC 60601-1-2:2014 Table 2 Requirements:

The GUS is intended for use in the electromagnetic environment specified below. The customer or the user of the GUS should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV Contact ±2 kV, ±4 kV, ±8 kV and ±15 kV Air	±8 kV Contact ±2 kV, ±4 kV, ±8 kV, ±15 kV Air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line(s) to line(s) ±0.5 kV, ±1 kV, ±2 kV line(s) to earth	±0.5 kV, ±1 kV line(s) to line(s) ±0.5 kV, ±1 kV, ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_T (100 % dip in U_T) for 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 70% U_T (30% dip in U_T) for 25 cycles 0% U_T (100% dip in U_T) for 5 seconds	<5% U_T (>95 % dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.
Power Frequency Magnetic Field (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Table 13. Table 2 Requirements—Electromagnetic Environment—Guidance

IEC 60601-1-2:2014 Table 4 Requirements:

The GUS is intended for use in the electromagnetic environment specified below.
The customer or the user of the GUS should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms for ISM bands.	3 Vrms	<p>Professional healthcare Environment</p> <p>WARNING: Portable and mobile RF communications equipment such as diathermy, electrocautery, and RFID equipment may affect the device. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the screen should be observed to verify normal operation while operating. The device may temporarily experience a disruption in function while near devices that emit strong radiated fields. Basic safety will not be impacted. If undesirable effects are observed, try one of the following:</p> <ol style="list-style-type: none"> 1. Reorient or relocate the other equipment. 2. The Uroflowmeter and Sensor are internally powered and will not be impacted by conducted emissions while running on battery. 3. While charging, connect the other equipment into an outlet on a circuit different from that to which the Uroflowmeter is connected to. 4. Consult Bright Uro for help.

<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p> <p>27 V/m 385 MHz</p> <p>28 V/m 450 MHz</p> <p>9 V/m 710/745/780 MHz</p> <p>28 V/m</p> <p>810/870/930 MHz</p> <p>28 V/m 1720/1845/1970 Mhz 28 V/m 2450 MHz</p> <p>9 V/m 5240/5500/ 5785 MHz</p>	<p>3 V/m</p>	<p>Professional Healthcare Environment</p> <p>WARNING: Portable and mobile RF communications equipment such as diathermy, electrocautery, and RFID equipment may affect the device. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the screen should be observed to verify normal operation while operating. The device may temporarily experience a disruption in function while near devices that emit strong radiated fields. Basic safety will not be impacted. If undesirable effects are observed, try one of the following:</p> <ol style="list-style-type: none"> 5. Reorient or relocate the other equipment. 6. The Uroflowmeter and Sensor are internally powered and will not be impacted by conducted emissions while running on battery. 7. While charging, connect the other equipment into an outlet on a circuit different from that to which the Uroflowmeter is connected to. 8. Consult Bright Uro for help.
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p>			
<p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the GUS is used exceeds the applicable RF compliance level above, the GUS should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the GUS</p> <p>b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 14. Table 4 Requirements—Electronic Environment—Guidance



WARNING: Per IEC 60601-1-2, ed 4.1, portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the Glean system Uroflowmeter or Sensor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.



NOTE: If the measured field strength at the location where the device is used exceeds the aforementioned compliance level, the device should be monitored to ensure it is functioning properly. Should unusual performance characteristics (i.e., no data acquisition when expected) be observed, additional measures may be required, such as changing the alignment or location of the device. Also, see the Troubleshooting section for help with troubleshooting.

Transmitter and Receiver Product Specifications:

The Uroflowmeter and Sensor devices intentionally transmit and receive RF energy for communication. The Uroflowmeter uses Wireless Power Transfer (WPT) for charging.	
Parameter	Specification
Tx/Rx Frequencies	2.402 to 2.480 GHz
Max Power Output	1.0 mW EIRP
Modulation	BLE1M-GFSK
WPT Frequency	100 kHz

Table 15. Transmitter and Receiver Product Specifications

Recommended separation distances between portable and mobile RF communications equipment and the Glean Urodynamics System			
The Glean Urodynamics System (GUS) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the GUS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the GUS as recommended below, according to the maximum output power of the communications equipment.			
			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz <i>dd = 1.17 pp</i>	80 MHz to 800 MHz <i>dd = 1.17 pp</i>	800 MHz to 2.5 GHz <i>dd = 2.33 pp</i>
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30
For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Table 16. Minimum Separation Distance

The clinical functions of the Glean UDS system are to wirelessly transmit pressure information (Sensor) and volume information converted from mass (Uroflowmeter) to achieve its intended use. The loss or degradation of this performance during normal use and a single fault condition does not result in unacceptable risk to the patient or operator and does not

affect the safety of the system. Therefore, Glean does not have essential performance and does not require recurrent testing to maintain basic safety.

APPENDIX E: END-USER SOFTWARE LICENSE AGREEMENT

EULA, Terms and Conditions: www.gleanuds.com/EULA

APPENDIX F: GLOSSARY

TERMS^{i, ii} USED IN GLEAN URODYNAMICS TESTING

NOTE: Notations in italics indicate usage specific to Bright Uro's Glean Urodynamics System protocols.

Abdominal leak point pressure: the intravesical pressure at which urine leakage occurs due to increased abdominal pressure in the absence of a detrusor contraction.

- *Measured in cmH₂O*

Acontractile detrusor: absence of detrusor contraction under Urodynamic evaluation.

Area under the curve: a calculation of the area contained by the curve of a urethral pressure profile.

Bladder pressure: (Pves, intravesical pressure) pressure within the bladder

- *Measured in cmH₂O*

Cystometry: the measurement of the pressure-volume relationship of the bladder during filling. Measurements obtained during cystometry include bladder sensations, compliance, bladder capacity and the presence or absence of detrusor overactivity (DO). The graphical recording of the bladder pressure and volume over time is referred to cystometrogram (CMG).

- *GUS bladder pressure is recorded as Pves in the Glean Web App*

Detrusor overactivity: characterized by involuntary detrusor contractions during the filling phase – either spontaneous or provoked.

Detrusor overactivity incontinence: incontinence due to an involuntary detrusor contraction.

Enuresis: involuntary loss of urine, usually subcategorized as nocturnal enuresis meaning involuntary loss of urine during sleep.

Filling phase: (storage phase) often used to describe the CMG portion of a Urodynamic examination, this phase ends upon to voiding.

First desire to void: during Urodynamics, the feeling that would lead the patient to pass urine at the next convenient moment, but voiding can be delayed.

- *This can be recorded as an event labeled "1st Desire" in the Glean Mobile App (Clinician)*

Frequency: the complaint of voiding too often by day.

Functional profile length: the length of the urethra along which the urethral pressure exceeds bladder pressure. o calculated within a UPP segment, displayed in mm, as length of continence zone

Hesitancy: difficulty initiating voiding.

Idiopathic detrusor overactivity: (formerly "detrusor instability") incontinence due to an involuntary detrusor contraction with no defined cause.

Incompetent urethral closure mechanism: when the urethra allows leakage of urine in the absence of a detrusor contraction.

Incontinence: the involuntary loss of urine. May be further defined as: stress incontinence, urge incontinence, mixed (both stress and urge) incontinence, nocturnal enuresis, and situational incontinence.

- *This can be recorded as an event labeled "Leak" in the Glean Mobile App (Clinician and Patient)*

International Continence Society: (ICS) "The primary interest of the International Continence Society is to study storage and voiding function of the lower urinary tract, its diagnosis and the management of lower urinary tract dysfunction, and to encourage research into pathophysiology, diagnostic techniques and treatment."

- *This group sets standards for Urodynamic testing that all Bright Uro training follows.*

Intravesical pressure: (Pves) pressure measured within the bladder. Note that pressure within the bladder can come from two sources – pressure from the abdomen (Pabd) and pressure from the muscle surrounding the bladder (Pdet). Formerly known as “total intravesical pressure”.

- *GUS intravesical pressure recorded in CMG, pressure/flow, and UPP tests, measured in cmH2O*

Intrinsic sphincter dysfunction: (ISD) is usually indicated by maximum urethral closure pressure less than 20 cmH2O pressure, or ALPP less than 60 cmH2O pressure.

Leak point pressure: (LPP, ALPP, VLPP, CLPP) the intravesical pressure at which involuntary urine leakage is noted during increased abdominal pressure, in the absence of a detrusor contraction.

- *Measured in cmH2O*

Lower urinary tract symptoms: (LUTS) these may include frequency, urgency, incontinence, nocturia, recurrent urinary tract infections, and many others.

Maximum cystometric capacity: (capacity) the volume at which the patient can no longer delay voiding. During Urodynamics, this is usually the point at which permission to void is given.

- *Measured in ml*

Maximum urethral pressure: (MUP) maximum pressure of the measured profile.

- *Measured in cmH2O*

Micturition study: a pressure/flow study. This study includes pressure measurements such as Pves and Uroflow measurements. This allows documentation of the relationship between the pressure generated during the voiding event and the resultant flow rate and pattern.

- *The results of the GUS pressure/flow study may be viewed in the Glean Web App*

Neuropathic detrusor overactivity: (formerly hyperreflexia) detrusor overactivity where there is a relevant neurological condition.

Nocturia: the complaint that patient has to wake from sleep during the night one or more times to void.

Nocturnal enuresis: the complaint of loss of urine during sleep.

Normal detrusor function: the detrusor allows the bladder to fill with little or no change in intravesical pressure, with no involuntary contractions despite provocation.

Permission to void: Time at which clinician allows patient to void as denoted by an annotation placed at time of reported sensation of bladder capacity, recommended by ICS to document when patient was told to allow voiding. This helps differentiate between contractions that are involuntary, and contractions that are voluntarily generated to initiate voiding.

- *This can be recorded as an event labeled “Voiding Event” in the Glean Mobile App (Clinician and Patient)*

Phasic detrusor overactivity: Intermittent detrusor overactivity which occurs during filling, which may or may not lead to incontinence.

Post-void residual: (PVR) the volume of urine left in the bladder after voiding.

- *This can be recorded as an event labeled “PVR” in the Glean Mobile App (Clinician)*

Retention: a non-painful bladder, which remains palpable or percussible after the patient has passed urine. Such patients may be incontinent.

Sensation: in Urodynamics, the reported sensations during testing such as first sensation, first desire, strong desire, and sense of reaching bladder capacity.

- *These can be recorded as events in the Glean Mobile App (Clinician)*

Stress urinary incontinence: (SUI) the symptom of a loss of urine associated with exertion, often with cough or sneeze. This is considered a complaint unless proven urodynamically, when it then is known as Urodynamic stress incontinence (formerly genuine stress incontinence).

Strong desire to void: described as the persistent desire to void without fear of leakage.

- *This can be recorded as an event labeled “Strong Desire” in the Glean Mobile App (Clinician).*

Terminal detrusor overactivity: a single involuntary detrusor contraction occurring at capacity, which cannot be suppressed and results in incontinence, usually resulting in emptying of bladder.

Uninhibited: acting without conscious inhibition – often used to describe a bladder contraction which the patient is unable to suppress.

Urethra: the tube leading from the bladder to the outside of the body.

Urethral pressure: (Pura) the pressure needed to just open a closed urethra.

- *Measured in cmH₂O*

Urethral pressure profile: (UPP) the pressures recorded throughout the length of the urethra, measured by withdrawing the Sensor at a slow known rate (recommended: 1mm/sec).

- *Measured in cmH₂O*

Urethral relaxation incontinence: leakage due to urethral relaxation in the absence of raised abdominal pressure or a detrusor contraction.

Urgency: the complaint of a sudden compelling desire to pass urine which is difficult to defer.

- *This can be recorded as an event labeled “Urge” in the Glean Mobile App (Clinician and Patient)*

Urge incontinence: symptom of incontinence associated with a strong compelling desire to void.

Urinary tract infection: finding of microbiological evidence of significant bacteriuria and pyuria usually accompanied by symptoms such as increased bladder sensation, urgency, frequency, dysuria, urgency urinary incontinence, and/or pain in the lower urinary tract.

Urodynamic stress incontinence: (formerly genuine stress incontinence, SUI, or stress incontinence) the involuntary leakage of urine during increased intravesical pressure, in the absence of a detrusor contraction.

Valsalva: the attempt to forcibly exhale with a closed glottis – often used to increase intra-abdominal pressure.

- *Used to provoke stress incontinence, can be recorded as an event labeled “Valsalva” in the Glean Mobile App (Clinician)*

Voiding phase: (emptying phase) often used to describe the portion of a Urodynamic evaluation that records both pressures and flow parameters during a voiding event, this would immediately follow the “filling phase”.

- *This can be recorded as an event labeled “Voiding Event” in the Glean Mobile App (Clinician and Patient)*

APPENDIX G: ACRONYMS

- CMG: Cystometrogram

- GUS: Glean Urodynamics System
- LUTD: Lower urinary tract dysfunction
- LUTS: Lower urinary tract symptoms
- MS: Micturition study
- PFS: Pressure/flow studies
- Pves: Intravesical pressure

- PVR: Post-void residual
- UF: Uroflowmetry
- UPP: Urethral pressure profile
- UTI: Urinary Tract Infection

13 REFERENCES

- i Abrams, Cardozo, Fall, Griffiths et al. The Standardisation of Terminology of Lower Urinary Tract Function: Report from the Standardisation Sub-committee of the International Continence Society. Neurourol Urodyn 21: 167–78, 2002.
- ii D’Ancona CD, Haylen BT, Oelke M, Herschorn S, Abranches-Monteiro L, Arnold EP, Goldman HB, Hamid R, Homma Y, Marcelissen T, Rademakers K, Schizas A, Singla A, Soto I, Tse V, de Wachter S. An International Continence Society (ICS) Report on the Terminology for Adult Male Lower Urinary Tract and Pelvic Floor Symptoms and Dysfunction. Neurourol Urodyn. 2019 DOI: 10.1002/nau.23897